



# Labculture® G4

### **Class II Biological Safety Cabinets**

The Most Certified Energy-efficient, Safe, and Ergonomic Biosafety Cabinet in the World





#### LABCULTURE® PLUS G4 (LA2-EU G4) CLASS II BIOSAFETY







#### **Airflow Sensor**

- Monitors real-time airflow for safety
- Alert the user if airflow is insufficient

ESCO

#### USB Port, RJ45, and Zero Volt Relay Contact

- USB Port / RJ45 to send operational information to Building Management System (BMS)
- Zero Volt Relay Contact to turn ON/OFF exhaust blower and signal the building alarm



### Centurion 7" Capacitive Touchscreen Controller

- Displays all safety information on one large screen
- Shows cabinet parameters with intuitive 3D illustration
- Easy to use menu, similar to Smart Phone Apps
- Large buttons, easy to operate when wearing gloves
- Self-guidance to users to deal with specific situations
- Centered and angled down for easy reach and viewing
- Optional: 21 CFR Part 11 Compliance



#### **Motorized Window** -

- Aerosol tight window seal
- Window automatically stops at safe operating height
- Conveniently move the window by fingertip

## Unique Stainless Steel and Glass Hybrid Wall (E-Series)

- Large corner radius for easy cleaning
- Easy to reach service fixture and outlets
- Stainless steel side wall is available (S-Series)
- Hole-free side glass for increased safety



#### **User-friendly Work Tray**

- Largest useable area in the market
- Recessed to contain spillage
- Sloped perimeter for easy cleaning
- Large, easy to clean tray handle



#### **Raised Arm Rest**

- Prevent grille blocking
- Comfortable working posture
- Durable stainless steel construction

Esco Labculture® Plus G4 Class II Biosafety Cabinet wailable in 3 feet, 4 feet, 5 feet, and 6 feet models.



#### **Ergonomic Work Zone**

- 10° angle to optimize user comfort, reduce glare, and maximize reach into the work area
- Brightly illuminated with >1200 lux (111 ft. cd)
- Industry-leading dimmable LED for optimum work comfort
- Airtight seal port for cable/tube exit protected by a negative pressure side wall

### **CABINET, FEATURING ADVANCED TOUCHSCREEN CONTROLLER**

#### **Dual Energy-efficient DC ECM Blower**

- The leading energy efficient Class II Biosafety Cabinet in the world with 70% energy savings compared to AC motor
- Stable airflow despite building voltage fluctuations and filter loading
- Standby mode to further reduce power consumption by 80%





Conventional AC Motor


Motor Voltage Fluctuation

Airflow Fluctuation

Esco DC Motor

Stable Motor Voltage

Stable Fluctuation



### **Advance Filtration System**

- 10x Filtration efficiency of HEPA filter
- Creates ISO Class 3 work zone instead of Industrystandard ISO Class 5
- Same 10 years filter life and replacement cost as HEPA filters



#### **Dimmable LED**

Save energy and optimize work comfort

#### **Removable Paper Catch**

- Prevent objects from being pulled into blower plenum
- Removable for easy cleaning
- Optional pre-filter can be fitted


#### **Tray Support Beams**

- Support work tray evenly for less vibration
- Cleaning holder to easily wipe the drain pan



#### **ISOCIDE™** Powder Coat

- Silver-ion impregnated powder coat
- Inhibits microbial growth to improve safety
- Prevents the plenum from becoming biohazard landfill



Certification								
Standards Compliance	Performance	Air Quality	Filtration	Electrical Safety				
	EN 12469, TUV Nord, Germany	ISO 14644.1, Class 3, Worldwide, US Fed Std 209E, Class 1 USA JIS 89920, Class 3, Japan	EN-1822 (H14), Europe IEST-RP-CC001, USA	EN-61010-1, Europe IEC61010-1, Worldwide				

#### Dynamic Chamber™

- Blower plenum and side walls are surrounded by negative pressure
- Prevent contaminants from escaping outside
- Positive PressureNegative Pressure




- H14-filtered air
- Unfiltered / potentially contaminated air
- Room air / Inflow air

#### **Cabinet Filtration System**


- Ambient air is pulled through front grille to create inflow, without going into the work surface. Inflow is joined by half of the downflow, to create front air curtain that is fine-tuned to create a large performance envelope. The combined air stream travels through the back air column towards the blower.
- Approximately ⅓ of the air in the common plenum is exhausted through the ULPA filter to the room. The remaining  $\frac{2}{3}$  of the air is passed through the downflow ULPA filter and into the work area as a vertical laminar flow air to create ISO Class 3 work surface and prevents cross contamination.
- Near the work surface, the downflow splits. About half goes to the front grille, and half goes to the rear grille. A small portion enters the the side capture zones to prevent dead air corners (small blue arrows).
- The design was optimized to give large performance envelope, that provides operator and product protection at wide Inflow and Downflow variation from the Nominal point.

#### **Performance Envelope Test of LA2-EU**



- Nominal Airflow
- Personnel / Product Protection
- Area of Personnel / Product Protection
- ▲ No Personnel /
  Product Protection
- Area of no Personnel / Product Protection

#### **Engineering Drawing**



- 1. Exhaust Airflow Sensor
- 2 LISB Port
- 3. Zero Volt Relay Contact
- 4. Power Inlet
- 5. Centurion Touchscreen Controller
- 6. Angled Sash Glass Window
- 7. Max. Sash Opening
- 8. Single-piece Stainless Steel Back Wall
- 9. Single-piece Stainless Steel Work Tray
- 10. Electrical Panel
- 11. Dimmable LED Lamp
- 12. IV Bar Retrofit Kit Provision
- 13. Electrical Outlet Provision
- 14. Service Fixture Retrofit Kit Provision
- 15. VHP Port In (Optional)
- 16. Tempered Side Glass
- 17. Stainless Steel Arm Rest
- 18. Exhaust Collar (optional)
- 19. Exhaust Filter
- 20. DC-ECM Blower
- 21. Downflow Filter
- 22. UV Lamp Provision
- 23. Cable Port (right side only)

		Options and A	ccessories				
Anti-blowback Valve 10 inches	304 Stainless Steel	ABBV-10S 5170354					
Exhaust Collar		ECO-F-LA2-3-G4 5171097	ECO-F-LA2-4-G4 5171098	ECO-F-LA2-5-G4 5171099	ECO-F-LA2-6-G4 5171100		
UV Lamp		UV-15A UV-30A 5170251 5170255					
IV Bar		IV-910 5170499	IV-1215 5170231	IV-1520 5170500	IV-1825 5170501		
Electrical Outlet	Direct Mounted	EO-HC 5170035					
	GFCI	EO-GFCI 5170071					
	EU SF-Gas-20 mm	SF-1G20 5170410					
	Solenoid Valve	SOLENOID GAS VALVE 230 VAC 5070086					
	EU SF-Vacuum-20 mm	SF-1V20 5170457					
Service Fixtures	EU SF-Air-20 mm	SF-1A20 5170502					
Service rixtures	EU SF-Nitrogen-20 mm	SF-1N20 5170503					
	EU SF-Water-20 mm	SF-1W20 5170458					
	US SF-Universal-20 mm	SF-2U22 5170504					
	Copper Piping for SF	CU-Pipe 5170026					
SPEU-A & SPEU-B	Support Stand 28" Height	SPEU-3A0 5131343	SPEU-4A0 5131344	SPEU-5A0 5131455	SPEU-6A0 5131456		
SPEU-A & SPEU-B	Support Stand 36" Height	SPEU-3B0 5131457	SPEU-4B0 5131458	SPEU-5B0 5131459	SPEU-6B0 5131460		
Pipette Storage Shelf		5260327					
Arm Rest Padding		MEWREST 5170127					
Foot Rest		FT-REST 5170492					
Laboratory Chair		ME-LD-AR360 1150006					
IQ OQ Protocol		9010179					







SF-1\_



UV-\_A-L



SF-2U\_



FT-REST





ME-LD-AR360





STA-\_













IQOQ



MEWREST

TECHNICAL SPECIFICATIONS							
Labculture® Class II 220-240 VAC, 50/60 Hz	Glass Side Walls		LA2-3E8-EU G4 2011682	LA2-4E8-EU G4 2011684	LA2-5E8-EU G4 2011686	LA2-6E8-EU G4 2011688	
	Stainless Steel Side	Walls	LA2-3S8-EU G4 2011683	LA2-458-EU G4 2011685	LA2-5S8-EU G4 2011687	LA2-6S8-EU G4 2011689	
Nominal Size		0.9 meter (3')	1.2 meter (4')	1.5 meter (5')	1.8 meter (6')		
External Dimensions* (W x D x H)	Without Optional B	ase Stand	1035 x 858 x 1570 mm (40.7" x 33.8" x 61.8")	1340 x 858 x 1570 mm (52.8" x 33.8" x 61.8")	1645 x 858 x 1570 mm (64.7" x 33.8" x 61.8")	1950 x 858 x 1570 mm (76.8" x 33.8" x 61.8")	
Internal Dimensions (W x D x H)		915 x 623 x 720 mm (36.0" x 24.5" x 28.3")	1220 x 623 x 720 mm (48.0" x 24.5" x 28.3")	1525 x 623 x 720 mm (60.0" x 24.5" x 28.3")	1830 x 623 x 720 mm (72.0" x 24.5" x 28.3")		
Usable Work Area		0.45 m² (4.8 sq.ft.)	Area: 0.62 m² (6.7 sq.ft.)	0.76 m² (8.2 sq.ft.)	0.93 m² (10.0 sq.ft.)		
Sash opening			200 mm (8")				
Maximum Sash Opening			470 mm (18.5")				
Average Airflow Velocity	Inflow		0.45 m/s (88.5 fpm)				
	Downflow		0.32 m/s (63 fpm)				
ULPA Filter Typical Efficiency			99.999% efficiency at 0.3 to 0.1 μm				
Sound Emission (dBA)*		54.8	55.0	57.5	58.0		
	Fluorescent Lamp Intensity		>1000 Lux (>93 foot-candles)				
Light Intensity	Optional UV Lamp		253.7 nm				
Electrical Requirements	Nominal Power Consumption			178 W	285 W	315 W	
	Heat Load			607 BTU/Hr	972 BTU/Hr	1075 BTU/Hr	
	Main body		Electro-galvanized steel with white oven-baked epoxy-polyester Isocide™ antimicrobial powder-coated finish, 1.5 mm (0.06") / 16 gauge thick				
Cabinat Construction	Work Zone		Stainless steel Type 304 with No.4 finish, 1.5 mm (0.06") / 16 gauge thick				
Cabinet Construction	Work Zone	E Series	UV absorbing tempered glass, 6 mm (0.2"), colorless and transparent				
		S Series	1.5 mm (0.06") 16 gauge stainless steel, type 304, with 4B finish				
Net Weight			300 Kg (661 lbs)	330 Kg (727.5 lbs)	370 Kg (816 lbs)	460 Kg (1014 lbs)	
Shipping Weight							
Shipping Dimensions, Maximum (W x D x H)							
Shipping Volume Dimensions (W x D x H)							

<sup>\*</sup>Electrical power consumption is an measurement of new unit with clean filter operated within nominal setpoint. Result per unit may vary.



LIFESCIENCES GROUP

## Improving Lives through Science

## ESCO.



- Animal Research Workstation
- Biological Safety Cabinet
  CO<sub>2</sub> Incubator
- Ducted Fume Hood
- Ductless Fume Hood
- Filtered Storage Cabinet
- Laboratory Centrifuge
- Laboratory Oven and Incubator
- Laboratory Refrigerator and Freezer
- Laboratory ShakerLaminar Flow Cabinet
- PCR Cabinet
- Powder Weighing Balance Enclosure
- Thermal Cycler
- Ultra-low Temperature Freezer



- Benchtop Incubator
- ART Workstation
- Anti-Vibration Table
- Gas Analyser











Pharmaceutical Manufacturing, Pharmacy Compounding, and Bioprocessing Tools



























**CRDMO Services** 

#### **ESCO LIFESCIENCES GROUP**

42 LOCATIONS IN 21 COUNTRIES ALL OVER THE WORLD



Follow us on social media, download our apps, and scan the QR code for more info.

























Esco Micro Pte. Ltd. • 21 Changi South Street 1 • Singapore 486 777 Tel +65 6542 0833 • mail@escolifesciences.com www.escolifesciences.com

Esco Technologies, Inc. • 903 Sheehy Drive, Suite F, Horsham, PA 19044, USA Tel: +1 215-441-9661 • eti.admin@escolifesciences.com

Esco Lifesciences Group Offices: Bangladesh | China | Denmark | Germany | Hong Kong | India | Indonesia | Italy | Japan | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam

