

Streptavidin-HRP Conjugate

Product Number: NEL750001EA

STORAGE CONDITIONS

Store at -20oC until first use. Store at 4oC after thawing. Do not re-freeze this product.

DESCRIPTION

Streptavidin-HRP (NEL750) is a highly sensitive detector conjugates. These conjugates bind to biotin labeled proteins in Western blot, ELISA and immunohistochemistry assays. They also bind to nucleic acids labeled with biotinylated nucleotide analogs in Southern and/or Northern blot, colony/plaque screen and in situ hybridization (ISH) assays.

When used in homologous DNA dot hybridizations with a biotin N6 labeled DNA probe, as little as 100 femtograms of nucleic acid are detectable on Kodak X-OMAT Blue Autoradiography film with Revvity Renaissance® Nucleic Acid Chemiluminescence Reagents in 2 hours or Revvity Renaissance CDP-Star® Nucleic Acid Chemiluminescence Reagents in 5-10 minutes.

These conjugates can also be used for in situ hybridization (ISH) or immunohistochemistry (IHC) assays as part of Tyramide Signal Amplification (TSA), whereby deposition of tyramide-linked detector molecules are catalyzed by horseradish peroxidase for subsequent amplified signal detection. Both the Streptavidin-HRP and Streptavidin-AP conjugates can be used, in conjunction with chromagenic substrates, as reporter molecules for biotin-tyramide deposition. Streptavidin-HRP can also be used as a catalyst for tyramide-linked molecule deposition in assays where nucleic acids (in situ) or proteins (IHC) are labeled with biotin.

SUGGESTED WORKING CONDITIONS

The following are suggested conjugate concentration ranges for various assays. Individual conditions may require some conjugate titration to obtain optimal signal clarity.

SA-HRP

Membrane assays (e.g., Southern, Northern, Western, colony/plaque): Dilute conjugate 1:12,500-1:25,000 in 1X Phosphate Buffered Saline (PBS), pH 7.1-7.4 (e.g., Dulbecco's PBS), 0.1% SDS, 0.05% Tween-20, blocking reagent (e.g., 0.5% Blocking Reagent) or equivalent.

Slide assays (i.e., ISH and IHC) including TSA: Dilute conjugate 1:1250-1:2500 in 0.10 M Tris-HCl, pH 7.5, 0.15 M NaCl, blocking reagent (e.g., 2% sheep serum) or equivalent.

STABILITY

Streptavidin-HRP and Streptavidin-AP are stable for at least 6 months from the date of shipment, when stored at 4°C

HAZARD WARNING

This product is considered to be non-hazardous. Although the product is classified as non-hazardous, we strongly recommend using prudent laboratory practices: Avoiding unnecessary contact, use of gloves, eye protection, lab coats, etc. when using this or any other laboratory reagent.

CDP-Star is a trade mark of Tropix, Inc. CDP-Star is protected under one or more of U.S. patents 5,326,882 and 4,931,569.

The information provided in this document is for reference purposes only and may not be all-inclusive. Revvity, Inc., its subsidiaries, and/or affiliates (collectively, "Revvity") do not assume liability for the accuracy or completeness of the information contained herein. Users should exercise caution when handling materials as they may present unknown hazards. Revvity shall not be liable for any damages or losses resulting from handling or contact with the product, as Revvity cannot control actual methods, volumes, or conditions of use. Users are responsible for ensuring the product's suitability for their specific application. REVVITY EXPRESSLY DISCLAIMS ALL WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, REGARDLESS OF WHETHER ORAL OR WRITTEN, EXPRESS OR IMPLIED, ALLEGEDLY ARISING FROM ANY USAGE OF ANY TRADE OR ANY COURSE OF DEALING, IN CONNECTION WITH THE USE OF INFORMATION CONTAINED HEREIN OR THE PRODUCT ITSELF

NEL750001EA-R REV01

The Revvity logo, consisting of the word "revvity" in a lowercase, sans-serif font.