

Axygen® LTS Compatible Tips



Axygen LTS Compatible Tips are ergonomically designed to reduce discomfort and optimized to work with Rainin LiteTouch pipettors, all with the quality you have come to expect from Corning. By reducing the force required to seal and eject, Axxygen LTS Compatible tips can help you avoid using excess force to obtain an effective seal between the pipettor and pipet tip. Axxygen LTS Compatible tips are available in your favorite packaging formats including Axxygen HybridRack™, a sustainable packaging choice. All LTS tips are composed of USP Class VI polypropylene and manufactured in an ISO 13485 and ISO 14001 United States based facility.

- ▶ RNase-/DNase-free
- ▶ Human gDNA-free
- ▶ PCR inhibition-free
- ▶ Nonpyrogenic
- ▶ USP Class VI polypropylene

20 µL LTS Compatible Extended Length, Graduated, Fine Point



Cat. No.	Packaging	Sterile	Quantity	Price
LTF-20-R-S	Racked, filtered	Yes	5x10x96 pcs	311 € + VAT
LTFER-20-R-S	HybridRack, filtered	Yes	5x10x96 pcs	311 € + VAT
LT-20-R-S	Racked	Yes	5x10x96 pcs	187 € + VAT
LT-20-R	Racked	No	5x10x96 pcs	162 € + VAT
LT-20-RL	Reload stack	No	4x10x96 pcs	100 € + VAT

200 µL LTS Compatible, Graduated, Fine Point



Cat. No.	Packaging	Sterile	Quantity	Price
LTF-200-R-S	Racked, filtered	Yes	5x10x96 pcs	311 € + VAT
LTFER-200-R-S	HybridRack, filtered	Yes	5x10x96 pcs	311 € + VAT
LT-200-R-S	Racked	Yes	5x10x96 pcs	187 € + VAT
LT-200-R	Racked	No	5x10x96 pcs	177 € + VAT
LT-200-RL	Reload stack	No	4x10x96 pcs	103 € + VAT
LT-200-TS	Tip station	No	5x5x192 pcs	196 € + VAT
LT-200-TS-S	Tip station	Yes	5x5x192 pcs	215 € + VAT

1000 µL LTS Compatible, Graduated, Fine Point



Cat. No.	Packaging	Sterile	Quantity	Price
LTF-1000-R-S	Racked, filtered	Yes	5x8x96 pcs	311 € + VAT
LTFER-1000-R-S	HybridRack, filtered	Yes	5x10x96 pcs	399 € + VAT
LT-1000-R-S	Racked	Yes	5x8x96 pcs	204 € + VAT
LT-1000-R	Racked	No	5x8x96 pcs	177 € + VAT
LT-1000-RL	Reload stack	No	4x8x96 pcs	100 € + VAT