



SARTORIUS ENTRIS BALANCES MID-YEAR CAMPAIGN

Entris® II essential line -30% off the price

BALANCE THAT FITS TO ANY LABORATORY - ENTRIS® II

IN SELECTION:

- Weighing capacity 60-8200 g
- Precision 0,1-1000 mg
- Internal calibration
- 12 programs
- RS232 connection
- GLP|GMP printout
- Verified and regular balances

FEATURES:

- Automatic calibration dependent on time and temperature (isoCAL)
- LED touch-screen
- USB-C connection
- Plug&Play fast and easy connection to PC, no extra software is needed
- BPT surface for high chemical resistance
- Special-modified glass draft shield to minimize weighing errors from static electricity



TAKE A CLOSER LOOK AT **ENTRIS® II** GENERATION ON SARTORIUS [PAGE](#) AND [BROCHURE](#).

EXAMPLES OF MODELS:

REGULAR MODELS

Article	Name/weighing capacity/precision	Regular price (€)*	Sale price (€)*
<u>BCE64i-1S</u>	Lab balance with int, Cal, 60 g, 0,1 mg	2180.00	1526.00
<u>BCE124i-1S</u>	Lab balance with int, Cal, 120 g; 0,1 mg	2310.00	1617.00
<u>BCE224i-1S</u>	Lab balance with int, Cal, 220 g; 0,1 mg	2420.00	1694.00
<u>BCE223i-1S</u>	Lab balance with int, Cal, 220 g; 1 mg	1450.00	1015.00
<u>BCE323i-1S</u>	Lab balance with int, Cal, 320 g; 1 mg	1590.00	1113.00
<u>BCE423i-1S</u>	Lab balance with int, Cal, 420 g; 1 mg	1700.00	1190.00
<u>BCE623i-1S</u>	Lab balance with int, Cal, 620 g; 1 mg	1940.00	1358.00
<u>BCE622i-1S</u>	Lab balance with int, Cal, 620 g; 10 mg	1320.00	924.00
<u>BCE822i-1S</u>	Lab balance with int, Cal, 820 g; 10 mg	1450.00	1015.00
<u>BCE1202i-1S</u>	Lab balance with int, Cal, 1200 g; 10 mg	1590.00	1113.00
<u>BCE2202i-1S</u>	Lab balance with int, Cal, 2200 g; 10 mg	1700.00	1190.00
<u>BCE3202i-1S</u>	Lab balance with int, Cal, 3200 g; 10 mg	1830.00	1281.00
<u>BCE4202i-1S</u>	Lab balance with int, Cal, 4200 g; 10 mg	1940.00	1358.00
<u>BCE6202i-1S</u>	Lab balance with int, Cal, 6200 g; 10 mg	2180.00	1526.00
<u>BCE2201i-1S</u>	Lab balance with int, Cal, 2200 g; 100 mg	1450.00	1015.00
<u>BCE5201i-1S</u>	Lab balance with int, Cal, 5200 g; 100 mg	1590.00	1113.00
<u>BCE8201i-1S</u>	Lab balance with int, Cal, 8200 g; 100 mg	1700.00	1190.00
<u>BCE6200i-1S</u>	Lab balance with int, Cal, 6200 g; 1 g	1450.00	1015.00
<u>BCE8200i-1S</u>	Lab balance with int, Cal, 8200 g; 1 g	1590.00	1113.00

*VAT will be added.

Campaign prices are valid until 31.07.2026.



VERIFIED MODELS

Article	Name/weighing capacity/precision	Regular price (€)*	Sale price (€)*
<u>BCE64i-1CEU</u>	Lab balance with int, Cal, 60 g, 0,1 mg	2680.00	1876.00
<u>BCE124i-1CEU</u>	Lab balance with int, Cal, 120 g; 0,1 mg	2810.00	1967.00
<u>BCE224i-1CEU</u>	Lab balance with int, Cal, 220 g; 0,1 mg	2920.00	2044.00
<u>BCE223i-1CEU</u>	Lab balance with int, Cal, 220 g; 1 mg	1940.00	1358.00
<u>BCE323i-1CEU</u>	Lab balance with int, Cal, 320 g; 1 mg	2060.00	1442.00
<u>BCE423i-1CEU</u>	Lab balance with int, Cal, 420 g; 1 mg	2180.00	1526.00
<u>BCE623i-1CEU</u>	Lab balance with int, Cal, 620 g; 1 mg	2420.00	1694.00
<u>BCE622i-1CEU</u>	Lab balance with int, Cal, 620 g; 10 mg	1830.00	1281.00
<u>BCE822i-1CEU</u>	Lab balance with int, Cal, 820 g; 10 mg	1940.00	1358.00
<u>BCE1202i-1CEU</u>	Lab balance with int, Cal, 1200 g; 10 mg	2060.00	1442.00
<u>BCE2202i-1CEU</u>	Lab balance with int, Cal, 2200 g; 10 mg	2180.00	1526.00
<u>BCE3202i-1CEU</u>	Lab balance with int, Cal, 3200 g; 10 mg	2310.00	1617.00
<u>BCE4202i-1CEU</u>	Lab balance with int, Cal, 4200 g; 10 mg	2420.00	1694.00
<u>BCE6202i-1CEU</u>	Lab balance with int, Cal, 6200 g; 10 mg	2680.00	1876.00
<u>BCE2201i-1CEU</u>	Lab balance with int, Cal, 2200 g; 100 mg	1940.00	1358.00
<u>BCE5201i-1CEU</u>	Lab balance with int, Cal, 5200 g; 100 mg	2060.00	1442.00
<u>BCE8201i-1CEU</u>	Lab balance with int, Cal, 8200 g; 100 mg	2180.00	1526.00
<u>BCE6200i-1CEU</u>	Lab balance with int, Cal, 6200 g; 1 g	1940.00	1358.00
<u>BCE8200i-1CEU</u>	Lab balance with int, Cal, 8200 g; 1 g	2060.00	1442.00

*VAT will be added.
Campaign prices are valid until 31.07.2026.

Ask more from our
Product managers:

nova@novabio.ee
nova@novabio.it



SARTORIUS

Authorized Partner

Sartorius and the Sartorius logo are registered trademarks of Sartorius AG or any of its affiliates.