



# SARTORIUS ENTRIS BALANCES MID-YEAR CAMPAIGN

**Entris® II essential line**  
**-30% off the price**

## BALANCE THAT FITS TO ANY LABORATORY - ENTRIS® II

### IN SELECTION:

- Weighing capacity 60-8200 g
- Precision 0,1-1000 mg
- Internal calibration
- 12 programs
- RS232 connection
- GLP|GMP printout
- Verified and regular balances

### FEATURES:

- Automatic calibration dependent on time and temperature (isoCAL)
- LED touch-screen
- USB-C connection
- Plug&Play fast and easy connection to PC, no extra software is needed
- BPT surface for high chemical resistance
- Special-modified glass draft shield to minimize weighing errors from static electricity



TAKE A CLOSER LOOK AT **ENTRIS® II** GENERATION ON SARTORIUS [PAGE](#) AND [BROCHURE](#).

## EXAMPLES OF MODELS:

### REGULAR MODELS

| Article                            | Name/weighing capacity/precision          | Regular price (€)* | Sale price (€)* |
|------------------------------------|---|--------------------|-----------------|
| <u><a href="#">BCE64i-1S</a></u>   | Lab balance with int, Cal, 60 g, 0,1 mg   | 2180.00            | 1526.00         |
| <u><a href="#">BCE124i-1S</a></u>  | Lab balance with int, Cal, 120 g; 0,1 mg  | 2310.00            | 1617.00         |
| <u><a href="#">BCE224i-1S</a></u>  | Lab balance with int, Cal, 220 g; 0,1 mg  | 2420.00            | 1694.00         |
| <u><a href="#">BCE223i-1S</a></u>  | Lab balance with int, Cal, 220 g; 1 mg    | 1450.00            | 1015.00         |
| <u><a href="#">BCE323i-1S</a></u>  | Lab balance with int, Cal, 320 g; 1 mg    | 1590.00            | 1113.00         |
| <u><a href="#">BCE423i-1S</a></u>  | Lab balance with int, Cal, 420 g; 1 mg    | 1700.00            | 1190.00         |
| <u><a href="#">BCE623i-1S</a></u>  | Lab balance with int, Cal, 620 g; 1 mg    | 1940.00            | 1358.00         |
| <u><a href="#">BCE622i-1S</a></u>  | Lab balance with int, Cal, 620 g; 10 mg   | 1320.00            | 924.00          |
| <u><a href="#">BCE822i-1S</a></u>  | Lab balance with int, Cal, 820 g; 10 mg   | 1450.00            | 1015.00         |
| <u><a href="#">BCE1202i-1S</a></u> | Lab balance with int, Cal, 1200 g; 10 mg  | 1590.00            | 1113.00         |
| <u><a href="#">BCE2202i-1S</a></u> | Lab balance with int, Cal, 2200 g; 10 mg  | 1700.00            | 1190.00         |
| <u><a href="#">BCE3202i-1S</a></u> | Lab balance with int, Cal, 3200 g; 10 mg  | 1830.00            | 1281.00         |
| <u><a href="#">BCE4202i-1S</a></u> | Lab balance with int, Cal, 4200 g; 10 mg  | 1940.00            | 1358.00         |
| <u><a href="#">BCE6202i-1S</a></u> | Lab balance with int, Cal, 6200 g; 10 mg  | 2180.00            | 1526.00         |
| <u><a href="#">BCE2201i-1S</a></u> | Lab balance with int, Cal, 2200 g; 100 mg | 1450.00            | 1015.00         |
| <u><a href="#">BCE5201i-1S</a></u> | Lab balance with int, Cal, 5200 g; 100 mg | 1590.00            | 1113.00         |
| <u><a href="#">BCE8201i-1S</a></u> | Lab balance with int, Cal, 8200 g; 100 mg | 1700.00            | 1190.00         |
| <u><a href="#">BCE6200i-1S</a></u> | Lab balance with int, Cal, 6200 g; 1 g    | 1450.00            | 1015.00         |
| <u><a href="#">BCE8200i-1S</a></u> | Lab balance with int, Cal, 8200 g; 1 g    | 1590.00            | 1113.00         |

\*VAT will be added.

Campaign prices are valid until 31.07.2026.



## VERIFIED MODELS

| Article                              | Name/weighing capacity/precision          | Regular price (€)* | Sale price (€)* |
|--------------------------------------|---|--------------------|-----------------|
| <a href="#"><u>BCE64i-1CEU</u></a>   | Lab balance with int, Cal, 60 g, 0,1 mg   | 2680.00            | 1876.00         |
| <a href="#"><u>BCE124i-1CEU</u></a>  | Lab balance with int, Cal, 120 g; 0,1 mg  | 2810.00            | 1967.00         |
| <a href="#"><u>BCE224i-1CEU</u></a>  | Lab balance with int, Cal, 220 g; 0,1 mg  | 2920.00            | 2044.00         |
| <a href="#"><u>BCE223i-1CEU</u></a>  | Lab balance with int, Cal, 220 g; 1 mg    | 1940.00            | 1358.00         |
| <a href="#"><u>BCE323i-1CEU</u></a>  | Lab balance with int, Cal, 320 g; 1 mg    | 2060.00            | 1442.00         |
| <a href="#"><u>BCE423i-1CEU</u></a>  | Lab balance with int, Cal, 420 g; 1 mg    | 2180.00            | 1526.00         |
| <a href="#"><u>BCE623i-1CEU</u></a>  | Lab balance with int, Cal, 620 g; 1 mg    | 2420.00            | 1694.00         |
| <a href="#"><u>BCE622i-1CEU</u></a>  | Lab balance with int, Cal, 620 g; 10 mg   | 1830.00            | 1281.00         |
| <a href="#"><u>BCE822i-1CEU</u></a>  | Lab balance with int, Cal, 820 g; 10 mg   | 1940.00            | 1358.00         |
| <a href="#"><u>BCE1202i-1CEU</u></a> | Lab balance with int, Cal, 1200 g; 10 mg  | 2060.00            | 1442.00         |
| <a href="#"><u>BCE2202i-1CEU</u></a> | Lab balance with int, Cal, 2200 g; 10 mg  | 2180.00            | 1526.00         |
| <a href="#"><u>BCE3202i-1CEU</u></a> | Lab balance with int, Cal, 3200 g; 10 mg  | 2310.00            | 1617.00         |
| <a href="#"><u>BCE4202i-1CEU</u></a> | Lab balance with int, Cal, 4200 g; 10 mg  | 2420.00            | 1694.00         |
| <a href="#"><u>BCE6202i-1CEU</u></a> | Lab balance with int, Cal, 6200 g; 10 mg  | 2680.00            | 1876.00         |
| <a href="#"><u>BCE2201i-1CEU</u></a> | Lab balance with int, Cal, 2200 g; 100 mg | 1940.00            | 1358.00         |
| <a href="#"><u>BCE5201i-1CEU</u></a> | Lab balance with int, Cal, 5200 g; 100 mg | 2060.00            | 1442.00         |
| <a href="#"><u>BCE8201i-1CEU</u></a> | Lab balance with int, Cal, 8200 g; 100 mg | 2180.00            | 1526.00         |
| <a href="#"><u>BCE6200i-1CEU</u></a> | Lab balance with int, Cal, 6200 g; 1 g    | 1940.00            | 1358.00         |
| <a href="#"><u>BCE8200i-1CEU</u></a> | Lab balance with int, Cal, 8200 g; 1 g    | 2060.00            | 1442.00         |

\*VAT will be added.  
Campaign prices are valid until 31.07.2026.

Ask more from our  
product manager:

**Heili Anni**

+372 53677740

[heili@novabio.ee](mailto:heili@novabio.ee)



**SARTORIUS**

Authorized Partner

Sartorius and the Sartorius logo are registered  
trademarks of Sartorius AG or any of its affiliates.