

Technical Datasheet

		EN ISO 13485 certified
Article - No:	06-695-5300	Revision: 5
Remarks		single use only

Specification

Product	Filter tip PP, premium surface, 100-1.250µl, 96 pcs/rack, transparent, moulded rings, np premium bio, sterile R, CE/IVD
Classification	CE/IVD
Dimensions	length: 98,20 mm (+/- 0,2 mm) inner diameter: 0,90 mm (+/- 0,2 mm)
Colour	transparent
Graduation	moulded rings
Nominal Capacity	100-1.250µl
Material Characteristics	polypropylene (PP), premium surface filter: high density polyethylene (PE)
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62 sterile R SAL 10-6
Purity Grade	free of PCR inhibitors free of human DNA free of RNase and DNase free of ATP free of pyrogens and endotoxins
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per inner cs: 576 pcs per carton: 3.456 pcs per pallet: 103.680
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C



Technical Datasheet

Date: 07.12.20

Page 2 of 2

	pal. not stackable
Approved by nerbe plus quality department	