

# Technical Datasheet

EN ISO 13485 certified		
<b>Article - No:</b>	<b>07-122-0073</b>	<b>Revision: 1</b>
<b>Remarks</b>	<b>Suitable for common pipettes e.g. Eppendorf, Biohit, Gilson or Brand</b>	<b>single use only</b>

## Specification

<b>Product</b>	Pipette tip PP "classic", 1-200µl, universal fit, yellow, np pcr ready, autoclavable up to 121°C, CE/IVD
<b>Classification</b>	CE/IVD
<b>Dimensions</b>	length: 52,60 mm (+/- 0,2 mm)
<b>Colour</b>	yellow
<b>Nominal Capacity</b>	1-200µl
<b>Material Characteristics</b>	polypropylene (PP)
<b>Temperature range</b>	-20°C - +90°C
<b>Autoclavability</b>	autocl. up to 121°C
<b>Product Characteristics</b>	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
<b>Purity Grade</b>	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
<b>Packaging</b>	pcs per bag: 1.000
	pcs per carton: 25.000
	pcs per pallet: 500.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable

**Approved by nerbe plus quality department**

00-000-0001/06.2020

