



Product Specification
Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,
LD-PE, transparent

Page 1



Product description

Order number	86.1171.010
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 155 x 15 mm, material: LD-PE, transparent, droplet size: 35 µl - 55 µl, injected graduation up to 3 ml, sterile, 10 piece(s)/peel-pack

Product characteristics

Graduation	yes
------------	-----

Size

Sample volume	3,5 ml
Width of product	15 mm
Length of product	155 mm
Drop size	35 µl - 55 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-07-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,
LD-PE, transparent

Page 2

Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product transparent

Purity & certification

Product in accordance with MDD or IVDD IVD product
CE certified CE
Purity standard sterile
Sterilisation Gas sterilisation
Lot or serial number yes

Packaging

Minimum order qty. 1200
Type of inner packaging peel-pack
Piece(s) / inner box 300
Piece(s) / case 1200
Piece(s) / pallet 33600
Depth of box 261 mm
Width of box 161 mm
Height of box 231 mm
Depth of case 533 mm
Width of case 333 mm
Height of case 250 mm
Case volume 0,0444 cbm
Weight of product 0,0015 kg
Weight of case 3,06 kg
EAN of inner box 4038917353144
EAN case 4038917353113

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-07-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT