

Technical Datasheet

| | | |
|----------------------|--------------------|------------------------|
| | | EN ISO 13485 certified |
| Article - No: | 07-962-0301 | Revision: 7 |
| Remarks | | single use only |

Specification

| | |
|--|--|
| Product | Gel-loading pipette tip PP, universal fit, 200µl, round tip, 96 pcs/rack, np pcr ready |
| Dimensions | length: 67,20 mm (+/- 0,2 mm) |
| Nominal Capacity | 200µl |
| Material Characteristics | polypropylene (PP) |
| Product Characteristics | biologically inert |
| | physiologically harmless |
| | no heavy metals according to EC/94/62 |
| | Food contact: raw material in compliance with regulation EC 1935/2004 |
| Purity Grade | free of PCR inhibitors |
| | free of human DNA |
| | free of RNase and DNase |
| | free of pyrogens |
| Packaging | pcs per inner cs: 768 |
| | pcs per carton: 3.072 |
| | pcs per pallet: 69.120 |
| Storage | dry condition: 40%-60% rel. humidity |
| | bags / cases orig. closed |
| | storage temp. +10°C to max. +30°C |
| | pal. not stackable |
| Approved by nerbe plus quality department | |