

# Technical Datasheet

		EN ISO 13485 certified
<b>Article - No:</b>	<b>04-072-0470</b>	<b>Revision: 4</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Reservoir PP, 120ml, low profile, 12 channels, transparent, np pcr ready
<b>Dimensions</b>	ANSI/SLAS 1
<b>Colour</b>	transparent
<b>Nominal Capacity</b>	120ml
<b>Material Characteristics</b>	polypropylene (PP)
<b>Temperature range</b>	-80°C - 121°C
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
<b>Purity Grade</b>	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
<b>Shelf - Life</b>	not limited at controlled storage conditions
<b>Packaging</b>	pcs per bag: 5
	pcs per inner cs: 10
	pcs per carton: 20
	pcs per pallet: 2.800
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable

Approved by nerbe plus quality department