

Technical Datasheet

Date: 22.06.20

Page 1 of 1

EN ISO 13485 certified		
Article - No:	05-003-0800	Revision: 3
Remarks	fits: 05-021-0000, 05-041-000 and 05-061-000	single use only

Specification

Product	Cap PE for analyzer cup, Ø16x8,5 mm, without cross slit, natural
Dimensions	diameter: 16,00 mm (+/- 0,2 mm)
	height: 8,50 mm (+/- 0,2 mm)
Colour	natural
Material Characteristics	polyethylene (PE)
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 1.000
	pcs per carton: 10.000
	pcs per pallet: 300.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	