

Technical Datasheet

EN ISO 13485 certified		
Article - No:	07-142-0000	Revision: 4
Remarks	suitable for: Eppendorf Reference, Socorex Calibra, Ratiopetta	single use only

Specification

Product	Pipette tip PP "classic", 500-2.500µl, transparent np pcr ready
Dimensions	length: 115,00 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	500-2.500µl
Material Characteristics	polypropylene (PP)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Packaging	pcs per bag: 250
	pcs per carton: 2.000
	pcs per pallet: 60.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

