

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-812-8071	Revision: 1
Remarks		single use only

Specification

Product	np GREEN SWAB flocked swab, very fine nylon head, ABS stick, breakpoint at 80mm, sterile R, CE2797
Classification	CE2797
Dimensions	Length: 150 mm (+/- 3mm) Width: 2,5 mm (+/- 0,2mm) Diameter head: 3 mm (+/- 1mm)
Material Characteristics	head: nylon stick: Acrylonitrile butadiene styren (ABS)
Product Characteristics	biologically inert physiologically harmless free of detectable phthalate acc. to REACH regulation
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per inner cs: 100 pcs per carton: 1.200 pcs per pallet: 129.600
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +2°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	