

# Technical Datasheet

EN ISO 13485 certified		
<b>Article - No:</b>	<b>11-162-0100</b>	<b>Revision: 1</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Mailing container PP, 60ml, Ø30x126 mm, with absorbing inlay, screw thread and print according to EN 829, without cap, natural, autoclavable up to 121°C
<b>Dimensions</b>	diameter: 30,00 mm (+/- 0,2 mm)
	height: 126,00 mm (+/- 0,2 mm)
<b>Colour</b>	natural
<b>Nominal Capacity</b>	60ml
<b>Material Characteristics</b>	polypropylene (PP)
<b>Temperature range</b>	-196°C - 121°C
<b>Autoclavability</b>	autocl. up to 121°C
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
<b>Packaging</b>	pcs per bag: 300
	pcs per carton: 300
	pcs per pallet: 7.200
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
<b>Approved by nerbe plus quality department</b>	