

Technical Datasheet

		EN ISO 13485 certified
Article - No:	11-163-0100	Revision: 1
Remarks		single use only

Specification

Product	Screw cap HDPE, for mailing container 30x126mm, with rim, natural
Colour	natural
Material Characteristics	high density polyethylene (HDPE)
Temperature range	-50°C – 70°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 300
	pcs per carton: 300
	pcs per pallet: 32.400
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	