

# Technical Datasheet

EN ISO 13485 certified		
<b>Article - No:</b>	<b>11-712-1010</b>	<b>Revision: 2</b>
<b>Remarks</b>	<b>suitable for formalin and chloroform</b>	<b>single use only</b>

## Specification

<b>Product</b>	Screw cap container PP, 1.250ml, Ø106x159 mm, with fitted screw cap PE, with EPE-pad, white, autoclavable up to 121°C (without pad)
<b>Colour</b>	white
<b>Nominal Capacity</b>	1.250ml
<b>Material Characteristics</b>	container: polypropylene (PP)
	cap: polypropylene (PP)
	pad: expanded polyethylene (EPE)
<b>Temperature range</b>	-80°C - 70°C
<b>Autoclavability</b>	autocl. up to 121°C (without pad)
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	flexible and soft
<b>Packaging</b>	pcs per cs: 35
	pcs per carton: 35
	pcs per pallet: 700
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
<b>Approved by nerbe plus quality department</b>	

00-000-0001/06.2020

