

Technical Datasheet

EN ISO 13485 certified	
Article - No:	02-351-6052
Remarks	Revision: 7 single use only

Specification

Product	Test tube PS with fitted screw cap PE, 12ml, Ø16x100 mm, round bottom, tube transparent, cap red, sterile R, max. RCF 3.000g
Dimensions	outer diameter: 16,00 mm (+/- 0,2 mm) height: 100,00 mm (+/- 0,2 mm)
Colour	tube transparent, cap red
Nominal Capacity	12ml
Material Characteristics	polystyrene (PS)
Temperature range	-10°C – 70°C
Centrifugation	max. RCF 3.000g
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months
Packaging	pcs per bag: 200 pcs per carton: 1.200 pcs per pallet: 28.800
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

