

Technical Datasheet

| | | |
|----------------------|--------------------|------------------------|
| | | EN ISO 13485 certified |
| Article - No: | 05-900-1040 | Revision: 2 |
| Remarks | | single use only |

Specification

| | |
|--|---|
| Product | Cryo vial PP, screw cap PE, 4ml, self-standing, transparent, grad., sterile R, autocl. |
| Dimensions | height: 73,3mm |
| Colour | transparent |
| Graduation | grad. |
| Nominal Capacity | 4ml |
| Material Characteristics | tube: polypropylene (PP) cap: polyethylene (PE) |
| Temperature range | -196°C - 121°C |
| Autoclavability | autocl. |
| Purity Grade | do not contain any pyrogens |
| Sterilisation Process | radiation sterilization SAL 10-6 |
| Packaging | pcs per bag: 100 pcs per carton: 500 pcs per pallet: 49.500 |
| Storage | dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +2°C to max. +30°C pal. not stackable |
| Approved by nerbe plus quality department | |