

Technical Datasheet

| | | |
|----------------------|--------------------|------------------------|
| | | EN ISO 13485 certified |
| Article - No: | 05-501-1020 | Revision: 3 |
| Remarks | | single use only |

Specification

| | |
|--|--|
| Product | Micro screw cap tube PP, skirted, 2ml, assembeled screw cap PP with o-ring, tube natural, cap transparent, sterile R, max. RCF 20.000g |
| Colour | tube natural, cap transparent |
| Nominal Capacity | 2ml |
| Material Characteristics | polypropylene (PP) |
| Temperature range | -196°C - 121°C |
| Centrifugation | max. RCF 20.000g |
| Purity Grade | free of RNase and DNase free of pyrogens |
| Sterilisation Process | radiation sterilization SAL 10-6 |
| Packaging | pcs per bag: 100 pcs per inner cs: 500 pcs per carton: 4.000 |
| Storage | dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +2°C to max. +30°C pal. not stackable |
| Approved by nerbe plus quality department | |

00-000-0001/06.2020

