

PureTrust™

ATP Intelligent Hygiene Monitoring System

PureTrust™ ATP Intelligent Hygiene Detection System is consisted by Intelligent Fluorescence Detector and Surface Swabs. Based on the principle of firefly luminescence, the system utilizes luciferase-luciferin's reaction system to quickly detect adenosine triphosphate (ATP) to monitor the hygiene status.

Product Features

Simple Operation and fast detection

Suitable for on-site testing.

Obtain accurate result within 30 seconds after sampling.

Easy for the staff to use with standardized operation.

Accurate and reliable

Test ATP presented in microbes and organic residues – the total amount of biological residues, to monitor cleaning effect of sanitation program.

Portable and flexible

Flexible swab design to achieve sampling even in tough area. The instrument can be operated with one hand.



Ordering information

Product Name	Catalog No.	Specification
PureTrust™ ATP Intelligent Fluorescence Detector	MF1000 premium	1 pc / box
PureTrust™ Surface Swab	PTS100	10 pcs / pack, 5 packs / box