

Technical Datasheet

		EN ISO 13485 certified
Article - No:	08-411-6000	Revision: 3
Remarks		single use only

Specification

Product	Inner container PS for sputum complete system, 24ml, Ø70x23mm, with fitted screw cap PE, container transparent, cap natural, sterile R
Dimensions	diameter: 70,00 mm (+/- 0,2 mm)
	height: 23,00 mm (+/- 0,2 mm)
Colour	container transparent, cap natural
Nominal Capacity	24ml
Material Characteristics	container: polystyrene (PS)
	cap: polyethylene (PE)
Temperature range	-20°C - 60°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Packaging	pcs per bag: 200
	pcs per carton: 200
	pcs per pallet: 9.600
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

